
anti-p-Stat3

Cat #: HM1356
Mouse monoclonal IgG
0.2 µg/µl, store at 4 °C

For research use only

BACKGROUND

Stats (signal transducers and activators of transcription) are transcription factors activated by peptide hormones, growth factors and cytokines, particularly interferons. STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators. STAT3 is activated through phosphorylation in response to various cytokines and growth factors including IFNs, EGF, IL5, IL6, HGF, LIF and BMP2. It mediates the expression of a variety of genes in response to cell stimuli, and thus plays a key role in many cellular processes such as cell growth and apoptosis.

SPECIFICITY

This antibody reacts specifically with Tyr-705 phosphorylated Stat3 of mouse, rat and human origin by Western blotting, immunoprecipitation, and immunohistochemistry

Molecular Weight of p-Stat3: 92 kDa.

IMMUNOGEN

A peptide corresponding to amino acid sequence containing phosphorylated Tyr-705 of human Stat3.

STORAGE

This antibody is stable for 12 months when stored at 2-8°C.

REFERENCES

1. Darnell JE Jr. 1997. STATs and gene regulation. *Science* 277:1630-5.
2. Darnell, J.E., Kerr, I.M., and Stark, G.R. 1994. Jak-STAT pathways and transcriptional activation in response to IFNs and other extracellular signaling proteins. *Science* 264: 1415-1421.
3. Zhong, Z., Wen, Z., and Darnell, J.E. 1994. Stat3: a STAT family member activated by tyrosine phosphorylation in response to epidermal growth factor and interleukin-6. *Science* 264: 95-98.
4. Hou, J., Schindler, U., Henzel, W.J., Ho, T.C., Brasseur, M., and McKnight, S.L. 1994. An interleukin-4-induced transcription factor: IL-4 Stat. *Science* 265: 1701-1706.
5. Pfeffer, L.M., Mullersman, J.E., Pfeffer, S.R., Murti, A., Shi, W. and Yang, C.H. 1997. STAT3 as an adapter to couple phosphatidylinositol 3-kinase to the IFNAR1 chain of the type I interferon receptor. *Science* 276: 1418-1420.
6. Nakashima, K., Yanagisawa, M., Arakawa, H., Kimura, N., Hisatsune, T., Kawabata, M., Miyazono, K. and Taga, T. 1999. Synergistic signaling in fetal brain by STAT3-Smad1 complex bridged by p300. *Science* 284, 479-482.
7. Bromberg, J.F., Wrzeszczynska, M.H., Devgan, G., Zhao, Y., Pestell, R.G., Albanese, C. and Darnell, J.E. Jr. 1999. Stat3 as an oncogene. *Cell* 98,295-303.

8. Zhang, T., Seow, K.T., Ong, C.T. and Cao, X. 2002. Interdomain interaction of Stat3 regulates its Src homology 2 domain-mediated receptor binding activity. *J. Biol. Chem.* 277, 17556-17563.
9. Leong, P.L., Andrews, G.A., Johnson, D.E., Dyer, K.F., Xi, S., Mai, J.C., Robbins, P.D., Gadiparthi, S., Burke, N.A., Watkins, S.F. and Grandis, J.R. 2003. Targeted inhibition of Stat3 with a decoy oligonucleotide abrogates head and neck cancer cell growth. *Proc. Natl. Acad. Sci. U.S.A.* 100, 4138-4143.

PRODUCTS FROM HYPROMATRIX, INC.**A. AntibodyArray™s:**

1. Signal Transduction AntibodyArray™
Catalog Number HM3000
2. Apoptosis AntibodyArray™
Catalog Number HM4000
3. Cell Cycle AntibodyArray™
Catalog Number HM5000

B. Staining AntibodyArray™s

1. Staining AntibodyArray™ I
Catalog Number HM8100
2. AntibodyArray Staining Apparatus
Catalog Number HM8000

C. Antibodies**1. HRP-conjugated antibodies**

- anti-phosphotyrosine
Catalog Number HM2040
- anti-phosphoserine
Catalog Number HM2070
- anti-phosphothreonine
Catalog Number HM2090

and more...

2. Primary antibodies

Hypromatrix offers a variety of high quality antibodies. For a complete list of antibodies and their specificities, please visit our web site at www.hypromatrix.com.

CONTACT

Hypromatrix, Inc.
100 Barber Avenue
Worcester, MA 01606
USA

Tel: 508-856-7900
Fax: 508-302-0748
Email: contact@hypromatrix.com
Web: www.hypromatrix.com